

Thermo Scientific Orion Star A111 pH Benchtop Meter

2013 Mid-year sales

from ~~RM 1,890.00~~
sale **RM 1,600.00** (save 15%)

Large, information screen clearly show pH or mV reading, temperature in °C or °F, meter mode, calibration information and power source

AUTO-READ™ function locks in the stable reading

Automatic recognition of USA / NIST and DIN buffers

Non-volatile memory holds up to 50 data points

Easily recall calibration slope data for procedure checks

3 year meter replacement warranty

The IP54-rated housing can handle splashes and can be wall-mounted

Included:

- Electrode
- Newly-designed probe holder
- Universal power adapter or
- Using AA batteries'

Up to 3 point pH calibration

Ex-Stock!

WORLD'S BEST SELLING BRAND

Come with:



**9157BNMD Triode 3-in-1
pH/ATC Electrode**



**60ml each of pH4.01, 7.00,
10.01 buffers and electrode
storage solution**



**Electrode arm with
redesigned holder**

NOTE: Promotion price above is based on Ex-Works.

Tel: 603-5621 1036

Fax: 603-5621 1037

Enquiry: sales@le.com.my

Website: www.le.com.my

**LABORATORY
EQUIPMENT**
SERVICING SCIENCE • PROVIDING SOLUTIONS

LE Laboratory Equipment Sdn Bhd